

IPR System

Solutions for most cases – The ContacEZ IPR Strip System provides the dental professional with accurate, easy-to-use devices to meet most patient IPR needs.

Contoured contacts - The high quality stainless steel strips are flexible to preserve and conform to the natural contours of the teeth during IPR. They are also tight enough to reach even the tightest contacts with ease.

Accurate – Paired with the Incremental Thickness Gauges, ContacEZ provides systematic, precise enamel removal without cutting lips, gums or cheeks.



The **ContacEZ IPR System** is an innovative set of dental instruments specifically crafted to make interproximal reduction safe, fast, and easy. These flexible dental strips are ergonomically designed, patient-friendly and each strip is marked with the thickness of the strip and whether it is single-sided (S) or double-sided (D). ContacEZ's interproximal reduction strips make IPR safer and more accurate without creating sharp corners or cutting lips and gums. Designed for use with clear aligners and orthodontic treatments, these strips prevent excessive enamel removal and maximize patient comfort and safety. Autoclavable and reusable.

IPR Strip System



Yellow Starter

0.07 mm
Single-Sided. Extra-Fine Diamond. Serrated
Prod No: 32108 8 pcs, 32132 32 pcs



Red Opener

0.12 mm
Double-Sided. Medium Diamond
Prod No: 32408 8 pcs, 32432 32 pcs



Blue Widener

0.15 mm
Double-Sided. Coarse Diamond
Prod No: 32508 8 pcs, 32532 32 pcs



Green Extra-Widener

0.20 mm
Double-Sided. Coarse Diamond
Prod No: 32608 8 pcs, 32632 32 pcs

IPR PLUS Strip System



Clear Opener

0.10 mm
Single-Sided. Med-Fine Diamond. Serrated
Prod No: 32208 8 pcs, 32232 32 pcs



Cyan Opener

0.15 mm
Single-Sided. Medium Coarse Diamond
Prod No: 32308 8 pcs, 32332 32 pcs



Purple Widener

0.25 mm
Double-Sided. Coarse Diamond
Prod No: 32708 8 pcs, 32732 32 pcs



Brown Widener

0.30 mm
Double-Sided. Coarse Diamond
Prod No: 32808 8 pcs, 32832 32 pcs



Mint Widener

0.35 mm
Double-Sided. Coarse Diamond
Prod No: 36608 8 pcs, 36632 32 pcs



Aqua Widener

0.40 mm
Double-Sided. Coarse Diamond
Prod No: 32708 8 pcs, 32732 32 pcs



Turquoise Widener

0.45 mm
Double-Sided. Coarse Diamond
Prod No: 36808 8 pcs, 36832 32 pcs



Teal Widener

0.50 mm
Double-Sided. Coarse Diamond
Prod No: 32908 8 pcs, 32932 32 pcs

Prod No Product Description

32008 IPR Assorted (0.07mm - 0.20mm) Small, 2 Yellow, 2 Red, 2 Dark Blue, 2 Green. 8 pcs

32024 IPR Rainbow Set (0.07mm - 0.50mm) 3 Yellow, 3 Red, 3 Blue, 3 Green, 2 Clear, 2 Cyan, 2 Purple, 2 Brown, 1 Mint, 1 Aqua, 1 Turquoise, 1 Teal. 24 pcs